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5. COMPONENTES

Hydrogen Peroxide, Purified Water, Disodium EDTA, Sodium Hydroxide.

6. EMBALAGEM

Kit - 12 frascos de 15ml

Concentrated solution for ultrasonic scaler and air prophy unit reservoir

INSTRUCTIONS FOR USE

Read all this leaflet carefully before you start using this product. Keep this leaflet, you may need to read it again.

1. WHAT IS ORALSAN OXYJET AND WHEN TO USE IT

ORALSAN OxyJet is concentrated 50% Hydrogen Peroxide and 1% Disodium EDTA to be used in ultrasonic scaler and air prophy unit reservoir. Each ORALSAN OxyJet vial is to be diluted in water to produce 500ml of 1.5% Hydrogen Peroxide solution to be used during scaling/air prophy procedures.

Such solution will treat drastically the intraoral environment reducing the bio-load of teeth, gums and aerosol developed during these procedures. ORALSAN OxyJet will facilitate scaling, stain removal and, at the same time, prevent the formation of calculus in equipment if not using deionized water.

2. HOW TO USE

Contains Hydrogen Peroxide 50%.

ATTENTION: wear gloves and protective goggles during solution preparation!

Eurogrin Pharma Ite
Filopimenos 3, Acharnes
13671- Greece
+30 210 2824437

Exclusive Distributor

IDS S.p.A.

Via S. Cristoforo 28/10
17100 Savona (SV) ITALY
Tel. +39 019 862080

www.ids dental.it

ORALSAN OxyJet must be diluted with deionized water before use. Use a proper plastic bottle with a capacity of at least 500 ml and add in this 485ml of deionized water. Open carefully a vial of ORALSAN OxyJet, remove the internal dropping tip of the vial and transfer the content into the plastic bottle. Mix well with a plastic spatula and transfer the solution into the tank of the ultrasonic device.

3. POSSIBLE SIDE EFFECTS

If you experience any of the following serious effects STOP using this product immediately: Urticaria, skin rash and itching, allergic reactions.

4. STORAGE

Keep out of the sight and reach of children. Do not store above 25°C. Store in original package. Do not use OxyWash MOUTHWASH after the expiry date stated on the label or carton.

5. INGREDIENTS

Hydrogen Peroxide, Purified Water, Disodium EDTA, Sodium Hydroxide.

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